



## SpecVette™ Price List - Jan 2008

Catalog Number	Unit Sale	Universal Adapter	Product Description	LIST PRICE
<b>SpecVette™ 250µm path length</b>				
SV250-25	25		SpecVette™ 250 µm pathlength, 5µL volume, 50 test kit	\$89.00
SV250-25A	25	Included	SpecVette™ 250 µm pathlength, 5µL volume, 50 test kit	\$129.00
<b>SpecVette™ 500µm path length</b>				
SV500-25	25		SpecVette™ 500 µm pathlength, 30µL volume, 50 test kit	\$99.00
SV500-25A	25	Included	SpecVette™ 500 µm pathlength, 30µL volume, 50 test kit	\$139.00
<b>SpecVette™ 1 mm path length</b>				
SV1000-25	25		SpecVette™ 1 mm pathlength, 60µL volume, 50 test kit	\$119.00
SV1000-25A	25	Included	SpecVette™ 1 mm pathlength, 60µL volume, 50 test kit	\$159.00
<b>Universal SpecVette™ Adapter</b>				
AD180-01	1		SpecVette™ (1cm x 1cm) holder, z-height 15mm One universal adapter accommodates any of the above SpecVettes	\$40.00

all prices in U.S. dollars, sales tax may apply, prices subject to change

z-height = 15mm

ALine can custom make SpecVettes™ to your specifications

**877-707-8575**  
[www.ALineInc.com](http://www.ALineInc.com)  
[info@ALineInc.com](mailto:info@ALineInc.com)